

# Development of an aberration corrected PEEM at the ALS

**Jun Feng**

**ESG, ALS**

**X-rays and magnetism workshop, Apr.14-18, 2003 Lake Tahoe**

# PEEM3 Project team



H.Padmore, A.MacDoeell,  
A.Scholl, M.Marcus,  
P.Schmid, A.Doran,  
T.Warwick,ESG

Weishi Wan, D.Robin,  
Accelerator Physics

R.Schlueter, R.Duarte,  
N.Kelez, K.Petermann,  
D.Munson, Engineering


S.Anders

E.Forest

D.H.Wei

Y.Wu




+ in the past, Mike Scheinfein



# Modeling the resolution




- Resolution evaluated by mapping secondary electron angle and energy distribution through fields
- electron distribution from sample in  $(\theta, \phi, E)$ 
  - Distribution in  $N(\theta) = \cos(\theta)$
  - Uniform distribution in  $\phi$
  - kinetic energy distribution  $N(E) = E/(E + w_f)^4$
- >Weighted “macro-particle”  
$$N(\theta, \phi, E) = \cos(\theta) \sin(\theta) E / (E + w_f)^4$$
- all order ray tracing technique used
- Detailed model check with SMART, PEEM2, charge ring +DA map method
- Modulation Transfer Function calculated by FFT method




# Aberration existed in electron lens just like in light optics lens



## Chromatic aberration



## Spherical aberration



$$r = \sqrt{(C_s \sin^3 \alpha)^2 + (C_c \delta \sin \alpha)^2 + (0.61\lambda / \sin \alpha)^2}$$


These aberrations limit PEEM resolution

# Electron mirror can compensate the aberration of electron lens




## Spherical aberrations



### Lens




### Mirror



## Chromatic aberrations



## Model checking with other theory




$$r = \theta \Delta f + \theta^3 C_3 + \theta^5 C_5 - \kappa \theta C_c - \kappa \theta^3 C_{3c} + \kappa^2 \theta C_{cc}$$

Aberration coefficients comparison

	$\Delta f$	$C_c$	$C_3$	$C_5$	$C_{3c}$	$C_{cc}$
Smart	$2.10^{-9}$	-9.8173m	-539,71m	-904640m	-3269.8m	-5.4317m
Our	$1.03.10^{-6}$	-9.829m	-539.84 m	-913725 m	-3268.5m	-5.387m
Diff	$\sim 0$	0.12%	0.02%	1%	0.04%	0.8%

# Model check with experiment: PEEM2 Resolution



**It begins with near unity modulation at low frequency and gradually becomes lower towards higher frequencies.**


**Contrast 1 means the image detail is perfectly maintained, while zero modulation transfer shows that it was completely lost, or not “seen” in the imaging process. Values between 0-1 indicate varying degrees of spatial details preserved.**

**50% resolution definition: 100nm**

**9% Rayleigh limit resolution: 38nm**

**MTF provides a continuum of unique ranking by which to judge a microscope’s resolution performance.**

# Model checking with experiment: PEEM2 Transmission









# Aberration correction




No mirror



With mirror



# Resolution for PEEM3




-All order ray trace results

-Diffraction limit included


-Both resolution and Transmission improved

# Aberration range




**Tetrode mirror  
covers the whole  
range for different  
operation conditions**

# PEEM3 layout




- high stability
- separator is in vacuum
- Two CCD



# Manipulator



# Dodecapole stigmator



# Electromagnetic dodecapole deflector







# 3D engineering modeling of Magnetic separator




Groove/ coil, water cooling  
through the wire






Electron trajectory



# PEEM3 beamline



# Zero order tracking for constant energy scale




$10^{-4}$ , stability  $<0.165\text{meV}$

# Floor Layout of MES /PEEM Sector 11.0

## Schedule

- microscope systems complete june 04
- commissioning with laser UV source
- EPU / beamline complete aug. 04




**Soft X-ray beamline with variable line space monochromator**

- 3 to 50 micron dia. Focus
  - bendable K-B mirrors
- dynamic measurement of photon energy
  - for high precision dichroism
  - zero order tracker


# PEEM2 - PEEM3 comparison



	PEEM2	PEEM3
Optics	Electrostatic lens	Electrostatic lens Electrostatic mirror Magnetic separator
Diagnostic	last image	movable pinhole, alignment PEEM
corrector	Octopole	Electromagnetic dodacapole
Resolution	20nm	5nm
Transmission	5%	>90%
@50nm		
Beamline	Bending 7.3.1.1	EPU 11.0.2
Relative	1	>1000
Flux density		



Schematic of PEEM3



Schematic of PEEM2